

**Tabel Usia Teknis Peralatan Medis(sumber AHA)**

No.	Medical Equipment	Life in Years.	Life in Units	Unit of Measure
1	Absorptiometer, Dual Photon, X-ray	8	200,800	exposure
2	Absorptiometer, SinglePhoton, X-ray	8	200,800	exposure
3	Aerator	12	24,096	hours
4	Air Compressor	15	30,120	hours
5	Air Sampler	10	20,080	hours
6	Alcohol Analyzer	8	16,064	hours
7	Amalgamator	6	12,048	hours
8	Amino Acid Analyzer	8	16,064	hours
9	Analyzer, Alcohol	8	16,064	hours
10	Analyzer, Amino Acid	8	16,064	hours
11	Analyzer, Antibody Susceptible	8	16,064	hours
12	Analyzer, Blod Gas/ph	8	16,064	hours
13	Analyzer, BUN	6	12,048	hours
14	Analyzer, Carbon Dioxide	5	10,040	hours
15	Analyzer, Carbon Monoxide	8	16,064	hours
16	Analyzer, Cell	10	20,080	hours
17	Analyzer, Clinical Chemistry	8	16,064	hours
18	Analyzer, Coagulation	5	10,040	hours
19	Analyzer, Defibrillator	8	16,064	hours
20	Analyzer, ECG Monitor	5	10,040	hours
21	Analyzer, Electrical Safety	8	16,064	hours
22	Analyzer, Electrolyte	7	14,056	hours
23	Analyzer, Electrosurgical	10	20,080	hours
24	Analyzer, Enzyme	5	10,040	hours
25	Analyzer, Esophageal Motility	7	14,056	hours
26	Analyzer, Glucose	7	14,056	hours
27	Analyzer, Hearing Aid	10	20,080	hours
28	Analyzer, Infusion Pump	5	10,040	hours
29	Analyzer, Lactate	7	14,056	hours
30	Analyzer, Lithium	8	16,064	hours
31	Analyzer, Middle Ear	10	20,080	hours
32	Analyzer, Nitrogen	8	16,064	hours
33	Analyzer, Nitrous Oxide	8	16,064	hours
34	Analyzer, Oxygen	10	20,080	hours
35	Analyzer, Pacemaker Function	8	16,064	hours
36	Analyzer, Platelet Aggregation	5	10,040	hours
37	Analyzer, Pulmonary Fuction	10	20,080	hours
38	Analyzer, Radioimmunoassay	8	16,064	hours
39	Analyzer, Sodium/Potassium	8	16,064	hours
40	Analyzer, Ultrasonic Unit	9	18,072	hours
41	Analyzer, Vision Function	10	20,080	hours
42	Anesthesia Apparatus	10	20,080	hours
43	Antibiotic Susceptibility Analyzer	8	16,064	hours
44	Apnea Monitor	7	61,320	hours
45	Arthroscope System	10	20,080	hours
46	Articulator, Dental	10	20,080	hours
47	Aspirator, Dental	10	20,080	hours
48	Aspirator, Surgical	10	20,080	hours
49	Aspirator, Thoracic	10	20,080	hours
50	Aspirator, Uterine	10	20,080	hours
51	Audiometer	10	20,080	hours
52	Audiometric Booth	15	131,400	hours
53	Auditory Function Screening Device	10	20,080	hours
54	Auscultoscope	10	20,080	hours
55	Autotransfusion Unit	8	16,064	hours

**Tabel Usia Teknis Peralatan Medis - Lanjutan**

<b>No.</b>	<b>Medical Equipment</b>	<b>Life in Years.</b>	<b>Life in Units</b>	<b>Unit of Measure</b>
56	Balance, Electronic	7	20,080	hours
57	Balance, Mechanical	10	20,080	hours
58	Basal Metabolism Unit	8	16,064	hours
59	Bassinet	12	105,120	hours
60	Bath, Paraffin, Clinical, Laboratory	10	20,080	hours
61	Bath, Paraffin, Physical Therapy	12	16,064	hours
62	Bath, Sitz	10	20,080	hours
63	Bath, Water, Electric	10	20,080	hours
64	Bath, Whirlpool	10	20,080	hours
65	Battery Charger	10	20,080	hours
66	Bed, Air Fluidized	10	20,080	hours
67	Bed, Birthing	12	105,120	hours
68	Bed, Circle, Electric	12	105,120	hours
69	Bed, Electric	12	105,120	hours
70	Bed, Hydraulic	12	105,120	hours
71	Bed, Pediatric	15	131,400	hours
72	Billirubinometer	8	16,064	hours
73	Biofeedback System	8	16,064	hours
74	Blast Cleaning Unit, Dental	10	20,080	hours
75	Blender, Clinical Laboratory	10	20,080	hours
76	Blood Cell Processor	10	20,080	hours
77	Blood Flow Detector, Ultrasonic	10	20,080	hours
78	Blood Gas/ph Analyzer	8	16,064	hours
79	Blood Glucose Monitor	8	70,080	hours
80	Blood Grouping System, Automated	9	18.072	hours
81	Blood Pressure Apparatus, Electronic	8	16,064	hours
82	Blood Recovery System	8	16,064	hours
83	Blood Warmer	8	16,064	hours
84	Booth, Audiometric	15	131.400	hours
85	Breast Pump	10	20,080	hours
86	Bronchoscope, Fiberoptic (FO)	6	12.048	hours
87	Bronchoscope, Rigid	10	20,080	hours
88	BUN Analyzer	6	12.048	hours
89	Cabinet, Dental (Sand Blasting )	10	20,080	hours
90	Cabinet, Warming	12	105.120	hours
91	Calibrator, Audiometer	8	16,064	hours
92	Calibrator-Analyzer	8	16,064	hours
93	Calibrator/Generator ECG	8	16,064	hours
94	Caloric Irrigation Unit	10	20,080	hours
95	Camera, Endoscopic	8	16,064	hours
96	Camera, Fundus	8	96,384	hours
97	Camera, GAMMA	8	80,320	hours
98	Camera, Identification (X-ray Film)	8	96,384	hours
99	Camera, Multi-Image	10	120,480	hours
100	camera, Radiographic Photospot	8	96,384	hours
101	Camera, Scintillation	8	96,384	hours
102	Camera, Surgical	8	321,280	hours
103	Camera, Video	8	16,064	hours
104	Capnograph	10	20,080	hours
105	Carbon Dioxide, Analyzer	5	10.040	hours
106	Carbon Dioxide, Monitor, Exhaled Gas	10	87,600	hours
107	Carbon Dioxide, Monitor, Transcutaneous	10	87,600	hours
108	Carbon Monoxide Analyzer	8	16,064	hours
109	Cardiac Output Unit	6	12.048	hours
110	Cart, Resuscitation	10	87,600	hours
111	Casting Unit, Dental	8	16.064	hours
112	Cell, Analyzer	10	20,080	hours
113	Cell, Saver	8	16,064	hours
114	Cell, Separator	8	16,064	hours
115	Cell, Washer	5	10.040	hours

No.	Medical Equipment	Life in Years.	Life in Units	Unit of Measure
116	Central Gas, Systems	10	87,600	hours
117	Centrifuge, Cell Washing	8	16,064	hours
118	Centrifuge, Laboratory, floor	10	20,080	hours
119	Centrifuge, Laboratory, High Speed	6	12,048	hours
120	Centrifuge, Refrigerated	8	16,064	hours
121	Centrifuge, Table Top	8	16,064	hours
122	Chair, Chamber, Hyperbaric	15	60,240	hours
123	Chair, Chloridometer	10	41,600	hours
124	Chair, Chromatography Equipment Gas	10	41,600	hours
125	Chair, Circulatory Assist Unit, Intra-Aortic Ballon	8	32,128	hours
126	Chair, Dental	10	41,600	hours
127	Chair, Examination/ Treatment	10	41,600	hours
128	Chair, Medical Instrument & Supply Sets	16	140,160	hours
129	Chair, Optometry	10	41,600	hours
130	Chair, Podiatric	10	41,600	hours
131	Charger, Battery	10	20,080	hours
132	Cleaner, Bedpan	10	20,080	hours
133	Cleaner, Dental, Steam	8	16,064	hours
134	Cleaner, Ultrasonic	10	20,080	hours
135	Cleaner, Vacuum	8	16,064	hours
136	Cleaner, Vacuum, Surgical	6	12,048	hours
137	Clinical Chemistry Analyzer	8	16,064	hours
138	Co-oximeter	10	20,080	hours
139	Coagulation Analyzer	5	10,040	hours
140	Coagulation Timer	8	16,064	hours
141	Coagulation, Sonic	10	20,080	hours
142	Collector, Dust	10	20,080	hours
143	Collimator, Radiographic	10	20,080	hours
144	Colonoscope	5	10,040	hours
145	Color Vision Tester	12	24,096	hours
146	Colposcope	7	14,056	hours
147	Compression Unit, Intermittent	10	20,080	hours
148	Compressor, Air	12	24,096	hours
149	Computerized Tomography (CT) Scanner prior to 1993	8	481,920	hours
150	CT Scanner, After 1993	5	502,000	hours
151	CT Scanner, After 1993 W/Upgrades	20	2,008,000	hours
152	Counter, Blood Cell	10	20,080	hours
153	Counter, Cell	5	10,040	hours
154	Counter, Colony	8	16,064	hours
155	Counter, GAMMA	8	16,064	hours
156	Counter, Particle	8	16,064	hours
157	Counter, Pill	10	20,080	hours
158	Counter, Platelet	8	16,064	hours
159	Counter, Scintillation	8	16,064	hours
160	Chruser, Syringe and Needle	10	20,080	hours
161	Cryogenic Blood Banking Equipment	10	20,080	hours
162	Cryostat	10	20,080	hours
163	Cryosurgical Unit	10	20,080	hours
164	Cryosurgical Unit, Ophthalmic	10	20,080	hours
165	Curring Unit, Dental, Acrylic	8	16,064	hours
166	Cutter, Bone	10	20,080	hours
167	Cutter, Cast	8	16,064	hours
168	Cutter-Aspirator, Vitrectomy	8	16,064	hours
169	Cutter-Vacuum, Orthopedic Cast	6	12,048	hours
170	Cystic-Fibrosis Screening Devices	10	20,080	hours
171	Cystometer	10	20,080	hours
172	Cystoscope	10	20,080	hours
173	Cystourethroscopy	10	20,080	hours
174	Darkroom, X-ray	6	12,048	hours
175	Defibrillator (w/proper battery replacement)	8	8,000	deep cycle charges
176	Defibrillator Analyzer	8	16,064	hours

No.	Medical Equipment	Life in Years.	Life in Units	Unit of Measure
177	Defibrillator / Monitor	8	16,064	hours
178	Demineralizer	8	16,064	hours
179	Densitometer, Recording	8	16,064	hours
180	Dental Articulator	10	20,080	hours
181	Dental Aspirator	10	20,080	hours
182	Dental Cabinet ( Sand Blasting )	10	20,080	hours
183	Dental Chair	10	20,080	hours
184	Dental Engine	8	16,064	hours
185	Dental Instrument Sterilizer	15	30,120	hours
186	Dental Operating and Treatment Unit	10	20,080	hours
187	Depth Perception Apparatus	10	20,080	hours
188	Dermatome	10	20,080	hours
189	Diathermy Unit	10	20,080	hours
190	Diathermy Unit, Ultrasonic	12	24,096	hours
191	Diluter	10	20,080	hours
192	Dispenser, Paraffin	10	20,080	hours
193	Distilling Apparatus	15	30,120	hours
194	Dosimeter, Radiation	10	20,080	hours
195	Dosimeter, Sound	10	20,080	hours
196	Drill, Bone	10	20,080	hours
197	Drill, Cranial	10	20,080	hours
198	Dryer, X-ray Film	10	20,080	hours
199	Dynamometer	15	30,120	hours
200	Dynamometer Exercise System	12	24,096	hours
201	ECG Monitor	10	87,600	hours
202	ECG Monitor Analyzer	5	10,040	hours
203	EEG Monitor	10	87,600	hours
204	Electrical Safety Analyzer	8	16,064	hours
205	Electro-oculograph	8	16,064	hours
206	Electrocardiograph, Interpretive	8	16,064	hours
207	Electrocardiograph, Multichannel	8	16,064	hours
208	Electrocardiograph, Three Channel	8	16,064	hours
209	Electrocautery Unit	10	20,080	hours
210	Electroencephalograph	10	20,080	hours
211	Electrolyte Analyzer	7	14,056	hours
212	Electromyograph	10	20,080	hours
213	Electronic Balance	7	14,056	hours
214	Electronystagmograph	8	16,064	hours
215	Electrophoresis Unit	8	16,064	hours
216	Electroretinograph	8	16,064	hours
217	Electrosurgical Analyzer	10	20,080	hours
218	Electrosurgical Unit	8	16,064	hours
219	Endoscope	10	20,080	hours
220	Endoscopic Camera	8	16,064	hours
221	Endoscopic Power Supply	10	20,080	hours
222	ENT Treatment Unit	12	48,192	hours
223	Enteral Feeding Pump	10	20,080	hours
224	Enzyme Analyzer	5	10,040	hours
225	Ergometer	10	20,080	hours
226	Esophageal Motility Analyzer	7	14,056	hours
227	Evoked Potential Unit	10	20,080	hours
228	Exerciser, Physical Therapy	15	30,120	hours
229	Eye Movement Monitor	10	20,080	hours
230	Fetal Heart Detector	8	70,080	hours
231	Fetal Monitor	10	87,600	hours
232	Fiberoptic Bronchoscope	6	12,048	hours
233	Fiberoptic Light Instruments	8	16,064	hours
234	Fibrometer	10	20,080	hours
235	Filter/Tank Unit, Pharmaceutical	6	12,048	hours
236	Fluorometer	10	20,080	hours

No.	Medical Equipment	Life in Years.	Life in Units	Unit of Measure
237	Freezer, Blood Cell	10	20,080	hours
238	Freezer, Blood Plasma	12	24,096	hours
239	Freezer, Laboratory	10	20,080	hours
240	Fundus Camera	8	96,384	exposures
241	Furnace, Dental	8	16,064	hours
242	GAMMA Camera	8	80,320	exposures
243	GAMMA Counter	8	16,064	hours
244	Gas Sampling Unit	10	20,080	hours
245	Gas Scavenger System, Anesthesia Unit	10	87,600	hours
246	Gastroscope, Rigid	5	10,040	hours
247	Gastroscope, Flexible	5	10,040	hours
248	Generator, Hydrogen	10	20,080	hours
249	Glucose Analyzer	7	14,056	hours
250	Grinding/Polishing Machine, Dental	8	16,064	hours
251	Headlight	12	24,096	hours
252	Hearing Aid Analyzer	10	20,080	hours
253	Heart Rate Monitor	10	87,600	hours
254	Heart-Lung Bypass system	8	70,080	hours
255	Heart Sealing Unit	10	20,080	hours
256	Heat Treating Unit, Ophthalmic	10	20,080	hours
257	Heater, Heat Treatment Pad	10	20,080	hours
258	Hemodialysis Unit	5	10,040	hours
259	Hemoglobinometer	10	20,080	hours
260	Hood, Chemical	10	60,240	hours
261	Hood, Dental Laboratory	10	60,240	hours
262	Hood, Fume	10	60,240	hours
263	Hood, Isolation, Laminar Air Flow	10	60,240	hours
264	Hood, Microbiological	12	72,288	hours
265	Humidifier	10	87,600	hours
266	Hydrocollator	10	20,080	hours
267	Hyperbaric Chamber	15	30,120	hours
268	Hypo/Hypothermia Unit	10	20,080	hours
269	Hypodermic Injection Unit	8	16,064	hours
270	Hypothermia Apparatus	10	20,080	hours
271	Identification Camera (X-ray Film)	8	96,384	exposures
272	Illuminator, X-ray, Motorized	8	16,064	hours
273	Immunodiffusion Equipment	10	20,080	hours
274	Incubator, Aerobic	10	87,600	hours
275	Incubator, Anaerobic	10	87,600	hours
276	Incubator, Infant	10	87,600	hours
277	Incubator, Infant, Transport	8	48,192	hours
278	Incubator, Test Tube, Laboratory	10	40,160	hours
279	Infusion Controller	10	20,080	hours
280	Infusion Pump	5	43,800	hours
281	Infusion Pump Analyzer	5	10,040	hours
282	Injector	10	20,080	hours
283	Injector, Angiographic	10	20,080	hours
284	Insufflator	10	20,080	hours
285	Irrigator	10	20,080	hours
286	Isokinetic Exercise Equipment	12	24,096	hours
287	Lactate Analyzer	7	14,056	hours
288	Lantern, Color Perception Testing	12	24,096	hours
289	Laparoscope	10	20,080	hours
290	Laser, Argon	10	20,080	hours
291	Laser, Carbon Dioxide	10	20,080	hours
292	Laser, Dye	10	20,080	hours
293	Laser, Krypton	10	20,080	hours
294	Laser, Nd:YAG	10	20,080	hours
295	Laser, Ophthalmic	10	20,080	hours
296	Lathe, Dental Laboratory	8	16,064	hours

No.	Medical Equipment	Life in Years.	Life in Units	Unit of Measure
297	Lavage Unit, Dental , Oral	8	16,064	hours
298	Lavege Unit, Surgical	8	16,064	hours
299	Lensometer	12	24,096	hours
300	Lift, Patient	10	20,080	hours
301	Ligth(Lamp), Bilirubin	10	87,600	hours
302	Ligth(Lamp), Dental	10	87,600	hours
303	Ligth(Lamp), Dental Activator	8	70,080	hours
304	Ligth(Lamp), Examination	10	87,600	hours
305	Ligth(Lamp), Fiberoptic	8	16,064	hours
306	Ligth(Lamp), Infrared	10	87,600	hours
307	Ligth(Lamp), Microscopic	10	40,160	hours
308	Ligth(Lamp), Phototherapy	10	87,600	hours
309	Ligth(Lamp), Surgical	10	87,600	hours
310	Ligth(Lamp), Ultraviolet	10	87,600	hours
311	Line Isolation Monitor	25	219,000	hours
312	Linear Accelerator	7	14,056	hours
313	Liquid Oxygen System	10	175,200	hours
314	Lithium Analyzer	8	16,064	hours
315	Lithotripter, Extracorporeal	5	30,120	electrodes
316	Lithotripter, Ultrasonic	5	30,120	electrodes
317	Loader, X-ray Film	6	12,048	hours
318	Magnet, Eye	10	20,080	hours
319	Magnetic Resonance Imaging 0,5 Tesla	8	1,124,480	slices
320	Magnetic Resonance Imaging 1,5 Tesla	8	1,124,480	slices
321	Magnetic Resonance Imaging Mobile	5	702,800	slices
322	Manikin	6	12,048	hours
323	Marker, Ocular	8	16,064	hours
324	Mechanical Balance	10	60,240	hours
325	Mercury Monitor/Detector	10	20,080	hours
326	Metal Detector, Magnetic	14	28,112	hours
327	Meter, Foot Candle	8	16,064	hours
328	Meter, Leakage,RF	8	16,064	hours
329	Microscope, Contact Lens	8	16,064	hours
330	Microscope, Electron	10	20,080	hours
331	Microscope, Laboratory	12	105,120	hours
332	Microscope, Operating	12	24,096	hours
333	Microscope, Phase Contrast	12	24,096	hours
334	Microscope, Stereo	12	105,120	hours
335	Microscope, Ultraviolet	12	24,096	hours
336	Microsurgical Instruments	10	20,080	hours
337	Microtome, Cryostat	10	20,080	hours
338	Microtome, Rotary	10	20,080	hours
339	Middle Ear Analyzer	10	20,080	hours
340	Mixer, Clinical Laboratory	10	20,080	hours
341	Mixer, X-ray Film Chemistry	10	20,080	hours
342	Mixer-Investor, Vacuum	8	16,064	hours
343	Moist Therapy Pack Conditioner	10	20,080	hours
344	Monitor, Apnea	7	61,320	hours
345	Monitor, Blood Glucose	8	70,080	hours
346	Monitor, Carbon Dioxide, Exhaled Gas	10	87,600	hours
347	Monitor, Carbon Dioxide, Transcutaneous	10	87,600	hours
348	Monitor, ECG	7	87,600	hours
349	Monitor, EEG	8	70,080	hours
350	Monitor, Eye Movement	8	70,080	hours
351	Monitor, Fetal Heart, Ultrasonic	7	61,320	hours
352	Monitor, Heart Rate	7	87,600	hours
353	Monitor, Oxygen, Transcutaneous	8	70,080	hours
354	Monitor, Pressure	8	70,080	hours
355	Monitor, Pulse	7	87,600	hours
356	Monitor, Radiation	8	70,080	hours

No.	Medical Equipment	Life in Years.	Life in Units	Unit of Measure
357	Monitor, Respiration	10	87,600	hours
358	Monitor, Temperature	8	70,080	hours
359	Monitor, Transcutaneous	7	61,320	hours
360	Nebulizer, Pneumatic	10	20,080	hours
361	Nebulizer, Ultrasonic	10	20,080	hours
362	Nephroscope	10	20,080	hours
363	Nitrogen Analyzer	8	16,064	hours
364	Nitrous Oxide Analyzer	8	16,064	hours
365	Oculoplethysmograph	8	16,064	hours
366	Ophthalmometer	11	22,088	hours
367	Organ Preservation System	10	20,080	hours
368	Osmometer	10	20,080	hours
369	Oven, Dental	8	16,064	hours
370	Oven, Drying	10	20,080	hours
371	Oximeter, Ear	10	20,080	hours
372	Oximeter, In Vitro	10	20,080	hours
373	Oximeter, Pulse	5	10,040	hours
374	Oxygen, Analyzer	10	20,080	hours
375	Oxygen, Blender	10	20,080	hours
376	Oxygen, Concentrator	8	16,064	hours
377	Oxygen, Meter	8	16,064	hours
378	Oxygenator	10	20,080	hours
379	Oxygen-Air Proportioner	10	20,080	hours
380	Pacemaker, Cardiac, External	5	10,040	hours
381	Pacemaker, Function Analyzer	8	43,800	hours
382	Package Sealer	10	20,080	hours
383	Pad, Circulating Fluid, Localized Heating/Cooling		20,080	hours
384	Paraffin Bath, Clinical, Laboratory	10	20,080	hours
385	Paraffin Bath, Physical Therapy	12	24,096	hours
386	PCA Pump	5	42,000	hours
387	Percussor	5	10,040	hours
388	Perimeter	10	20,080	hours
389	Peritoneal Dialysis Unit	8	16,064	hours
390	ph Meter	10	20,080	hours
391	Phonocardiograph	8	16,064	hours
392	Phorometer	12	24,096	hours
393	Phoropter	12	24,096	hours
394	Photometer, Flame	8	16,064	hours
395	Photometer, Reflectance	10	20,080	hours
396	Phototachometer	5	10,040	hours
397	Phototherapy Unit	10	20,080	hours
398	Phototherapy Unit, Hyperbilirubinemia	8	16,064	hours
399	Physiologic Monitor System, Acute Care	7	61,320	hours
400	Physiologic Monitor System, Cardiac Catheterization	7	61,320	hours
401	Physiologic Monitor System, Neonatal	7	61,320	hours
402	Physiologic Monitor System, Strees Exercise	7	61,320	hours
403	Pill Counter, Pharmaceutical	10	40,160	hours
404	Pipette, Automatic	10	20,080	hours
405	Platelet Aggregation Analyzer	5	10,040	hours
406	Plethysmograph	10	20,080	hours
407	Pneumatic Tester	8	16,064	hours
408	Pneumothorax Unit	12	24,096	hours
409	Porcelain Furnace, Dental	10	20,080	hours
410	Pressure Monitor	10	87,600	hours
411	Pressure Pad, Alternating	10	20,080	hours
412	Processor, X-ray Film, Automatic, Dental	8	32,128	hours
413	Processor, X-ray Film, Automatic, Full Size	6	24,096	hours
414	Processor, X-ray Film, Automatic, Table Top	7	28,112	hours
415	Projector, Chart, Eye	10	40,160	hours
416	Projector, Cine	10	20,080	hours

No.	Medical Equipment	Life in Years.	Life in Units	Unit of Measure
417	Projector, Microscope	12	48,192	hours
418	Projector, Ophthalmic	10	40,160	hours
419	Prophylactic Unit, Dental, Ultrasonic	6	12,048	hours
420	Pulmonary Function Analyzer	10	20,080	hours
421	Pulse Monitor	10	87,600	hours
422	Pulse Oximeter	5	10,040	hours
423	Pump, Breast	10	20,080	hours
424	Pump, Enternal Feeding	10	20,080	hours
425	Pump, Infusion	5	43,800	hours
426	Pump, PCA	5	42,000	hours
427	Pump, Suction	10	20,080	hours
428	Pump, Suction, Oxytocin	10	20,080	hours
429	Pump, Vacuum	10	20,080	hours
430	Pump, Volumetric	10	20,080	hours
431	Radiation Counter	8	16,064	hours
432	Radiation Monitor	8	16,064	hours
433	Radiographic Photospot Camera	8	96,384	exposure
434	Radiographic Table Systems	10	240,960	overhead tube exp.
435	Radiographic Unit, Chest	8	192,768	overhead tube exp.
436	Radiographic Unit, Dental	8	192,768	overhead tube exp.
437	Radiographic Unit, General Purpose under 500 mA	12	289,152	overhead tube exp.
438	Radiographic Unit, General Purpose over 500 mA	8	192,768	overhead tube exp.
439	Radiographic Unit, Mammographic	10	240,960	overhead tube exp.
440	Radiographic Unit, Mobile	8	192,768	overhead tube exp.
441	Radiographic Unit, Therapeutic	8	192,768	overhead tube exp.
442	Radiographic Unit, Tomographic Computer Assisted	8	481,920	slices
443	Radiographic Unit, Tomographic Not Computer Assisted	10	240,960	overhead tube exp.
444	Radiographic/Fluoroscopic Unit, General Purpose	10	240,960	overhead tube exp.
445	Radiographic/Fluoroscopic Unit, Mobile	8	192,768	overhead tube exp.
446	Radioimmunoassay Analyzer	8	16,064	hours
447	Radiometer	8	16,064	hours
448	Radiometer, Ultraviolet	8	16,064	hours
449	Radiotherapy Unit	8	16,064	hours
450	Recorder, Chart	10	20,080	hours
451	Recorder, Echocardiographic	8	16,064	hours
452	Recorder, Long Term	6	12,048	hours
453	Recorder, Long Term ECG, Portable	5	10,040	hours
454	Recorder, Tape	8	16,064	hours
455	Recorder, Video Disk	8	16,064	hours
456	Refractometer	10	20,080	hours
457	Refractor, Ophthalmologic	10	20,080	hours
458	Refrigerator, Cold Pack Therapy	12	105,120	hours
459	Refrigerator, Biologic	12	105,120	hours
460	Refrigerator, Blood	12	105,120	hours
461	Refrigerator, Blood Bank	10	87,600	hours
462	Refrigerator, Laboratory	12	105,120	hours
463	Refrigerator, Morgue	12	105,120	hours
464	Refrigerator/Freezer	12	105,120	hours
465	Regulator, High Pressure, Gas	10	87,600	hours
466	Renal Transport Unit	8	16,064	hours
467	Respiration Monitor	10	20,080	hours
468	Resuscitator, Pulmonary, Gas	6	12,048	hours
469	Resuscitator, Pulmonary, Manual	10	12,048	hours
470	Rotator	10	20,080	hours
471	Sample, Air	10	20,080	hours
472	Sanitizer	12	48,192	hours
473	Saw, Autopsy	10	105,120	hours
474	Saw, Bone	10	20,080	hours
475	Saw, Surgical	10	20,080	hours
476	Scale, Autopsy	12	105,120	hours



No.	Medical Equipment	Life in Years.	Life in Units	Unit of Measure
477	Scale, Bed	10	20,080	hours
478	Scale, Floor	14	28,112	hours
479	Scale, Infant	15	30,120	hours
480	Scale, laboratory	8	16,064	hours
481	Scanner, Computerized Tomography (CT) prior to 1993	8	481,920	slices
482	Scanner, CT, After 1993	5	502,000	slices
483	Scanner, CT, After 1993 W/Upgrades	20	2,008,000	slices
484	Scanner, Long-term Recording, ECG	7	14,057	hours
485	Scanner, Ultrasonic, General Purpose	8	16,064	hours
486	Scintillation Camera	8	96,384	exposures
487	Scintillation Counter	8	16,064	hours
488	Sealing Machine	10	20,080	hours
489	Sealing Machine, Unit Dose	10	20,080	hours
490	Sensitometer, Radiographic	10	20,080	hours
491	Shaker, Laboratory	10	20,080	hours
492	Shapener, Dental	7	14,056	hours
493	Shapener, Microtome Knife	10	20,080	hours
494	Simulator, Arrhythmia	7	14,056	hours
495	Simulator, ECG	10	20,080	hours
496	Simulator, Medical Functions	8	16,064	hours
497	Simulator, Radiotherapy	10	20,080	hours
498	Sink, Portable	12	72,288	hours
499	Sink, Surgical Scrub	12	140,160	hours
500	Sitz Bath	10	20,080	hours
501	Slide Stainer	10	20,080	hours
502	Slit Lamp	10	20,080	hours
503	Sodium/Potassium Analyzer	8	16,064	hours
504	Spectrometer, Mass	8	16,064	hours
505	Spectrophotometer	8	16,064	hours
506	Sphygmomanometer, Electronic	8	16,064	hours
507	Spirometer	10	20,080	hours
508	Stapler, Automatic, Suture	5	10,040	hours
509	Sterilizer, Agar	10	40,160	hours
510	Sterilizer, Dental Instrument	15	30,120	hours
511	Sterilizer, Ethylene Oxide	16	32,128	hours
512	Sterilizer, Inoculating Loop	15	30,120	hours
513	Sterilizer, Steam (16X16X26)	12	36,000	hours
514	Sterilizer, Steam (20X20X38)	12	36,000	hours
515	Sterilizer, Steam (24XX Series)	15	45,000	hours
516	Sterilizer, Steam (26XX Series)	15	45,000	hours
517	Sterilizer, Steam, Laboratory	12	9,000	hours
518	Sterilizer, Steam, Table Top	10	30,000	hours
519	Sterilizer, Steam/Ethylene Oxide	16	32,128	hours
520	Stethoscope, Electronic	10	20,080	hours
521	Stethoscope, Ultrasonic	7	14,056	hours
522	Stimulator, Bone Growth	10	20,080	hours
523	Stimulator, Nerve	10	20,080	hours
524	Stimulator, Visual	18	36,144	hours
525	Stirrer	10	20,080	hours
526	Stretcher	15	131,400	hours
527	Suction Apparatus	10	20,080	hours
528	Suction Pump	10	20,080	hours
529	Suction/Pressure Unit	10	20,080	hours
530	Surgical Aspirator	10	20,080	hours
531	Surgical Camera	8	321,280	exposures
532	Syringe Pump, Oxytocin	10	20,080	hours
533	Table, Examination/Treatment	15	131,400	hours
534	Table, Fracture	15	131,400	hours
535	Table, Nuclear Medicine	10	40,160	hours
536	Table, Obstretical	20	175,200	hours

No.	Medical Equipment	Life in Years.	Life in Units	Unit of Measure
537	Table, Operating	15	131,400	hours
538	Table, Orthopedic	10	87,600	hours
539	Table, Physical Therapy	15	131,400	hours
540	Table, Urological	15	131,400	hours
541	Tachistoscope	10	20,080	hours
542	Tachometer, Strobe	8	16,064	hours
543	Telemetry Unit	10	20,080	hours
544	Temperature Monitor	8	70,080	hours
545	Tent, Oxygen	8	70,080	hours
546	Test Set, Electrosurgical	8	16,064	hours
547	Tester, Color Vision	12	24,096	hours
548	Tester, Current Leakage	8	16,064	hours
549	Tester, Defibrillator	8	16,064	hours
550	Tester, Pulp	8	16,064	hours
551	Thermometer, Electronic	5	10,040	hours
552	Thermoregulator	8	16,064	hours
553	Thoracic Aspirator	10	20,080	hours
554	Thyroid Uptake Systems	10	20,080	hours
555	Timer	10	20,080	hours
556	Timer, Coagulation	8	16,064	hours
557	Tissue Embedding Equipment	10	20,080	hours
558	Tissue Processor	10	20,080	hours
559	Titration	10	20,080	hours
560	Tonometer	10	20,080	hours
561	Tourniquet, Pneumatic	10	20,080	hours
562	Traction Unit	10	20,080	hours
563	Transcutaneous Monitor	7	87,600	hours
564	Transilluminator	12	61,320	hours
565	Treadmill	8	16,064	hours
566	Turning Frame, Orthopedic	10	20,080	hours
567	Tympanometer	10	20,080	hours
568	Ultrasonic Blood Flow Detector	10	20,080	hours
569	Ultrasonic Unit, Analyzer	9	18,072	hours
570	Ultrasonic Unit, Dental Prophylaxis	6	12,048	hours
571	Ultrasonic Unit, Diagnostic	12	24,096	hours
572	Ultrasonic Unit, Therapeutic	7	14,056	hours
573	Unit Dose Packing System	10	20,080	hours
574	Urinary Collection Unit, Precision Measurement	8	16,064	hours
575	Urodynamic Measurement System	7	14,056	hours
576	Uterine Aspirator	10	20,080	hours
577	Vacuum Mixing Device, Dental	8	16,064	hours
578	Vacuum Pump	10	20,080	hours
579	Vacuum System, Dental	12	24,096	hours
580	Vaporizer, Anesthesia Unit	10	20,080	hours
581	Vectorcardiograph	7	14,056	hours
582	Ventilator, Air Evacuation	10	10,000	hours
583	Ventilator, Anesthesia Unit	10	10,000	hours
584	Ventilator, Neonatal	10	10,000	hours
585	Ventilator, Volume	10	10,000	hours
586	Vibrator, Mixer, Dental	10	20,080	hours
587	Videa, Camera	8	16,064	hours
588	View Box, Rh Typing	10	87,600	hours
589	View Box, Vision Function Analyzer	10	87,600	hours
590	View Box, X-ray	12	105,120	hours
591	View Box, X-ray, Motorized	10	20,080	hours
592	Vision Function Analyzer	10	20,080	hours
593	Vision Tester, Color	12	24,096	hours
594	Vitrectomy Unit	7	14,056	hours
595	Volumetric Pump	10	20,080	hours
596	Warming Cabinet	12	72,288	hours

No.	Medical Equipment	Life in Years.	Life in Units	Unit of Measure
597	Warmer Radiant, Infant	10	20,080	hours
598	Warmer, Blood	8	16,064	hours
599	Warmer, Slide, Histology	10	20,080	hours
600	Warmer, Solution	12	24,096	hours
601	Washer, Bedpan	10	40,160	hours
602	Washer, Labware	8	32,128	hours
603	Washer, Surgical Instruments	8	32,128	hours
604	Washer/Sterilizer	8	48,000	hours
605	Water Bath, Electric	10	40,160	hours
606	Welding Unit, Orthodontic	8	16,064	hours
607	Wheelchair, Manual	6	52,560	hours
608	Wheelchair, Powered (with proper battery maintenance)	5	43,800	hours
609	Whirlpool Bath	10	20,080	hours
610	X-ray Absorptiometer, Single Photon	8	200,800	exposures
611	X-ray Absorptiometer, Dual Photon	8	200,800	exposures
612	X-ray Calibration/Verification System	8	16,064	hours
613	X-ray Chemical Mixer	10	20,080	hours
614	X-ray Film Changer, Cassette	10	20,080	hours
615	X-ray Film Dryer	8	16,064	hours
616	X-ray Film Duplicator	8	16,064	hours
617	X-ray Film Headling Equipment, Automatic, Daylight	6	12,048	hours
618	X-ray Film Identification Camera	8	96,384	exposures
619	X-ray Film Processor	8	16,064	hours
620	X-ray Film Processor, Automatic, Cline	6	12,048	hours
621	X-ray Film Processor, Automatic, Dental	8	16,064	hours
622	X-ray Film Processor, Automatic, Full Size	6	12,048	hours
623	X-ray Film Processor, Automatic, Table Top	7	14,056	hours

